



Kidney Biopsy Instruction Form

Section	Contents	Page
1	Contents of kidney biopsy kit	page 1
2	Tissue handling	page 2
3	Packing and return shipping instructions	page 4
4	Reporting of results and billing	page 5

1. Contents of the kidney biopsy kit

Please check all contents upon receipt. If any contents are missing or damaged, please contact:
MyLav Laboratory (info@mylav.net)

- a. Forms/instructions
 - Kidney biopsy instruction form (this form)
- b. Fixatives
 - Fixative for electron microscopy (TEM): Glutaraldehyde
 - Fixative for immunofluorescence (IF): Michel's transport medium
 - Fixative for light microscopy (LM): 4% buffered formalin



2. Tissue handling

Key aspects:

- Keep fixative for TEM (Glutaraldehyde) and IF (Michel's Medium) cool (+4° Celsius, not frozen).
- Keep tissue moist using physiologic saline and process biopsies as soon as possible (within 5 minutes) after collection/excision.
- During the process of partitioning and fixing of the biopsies, do not risk cross-contaminating tissue samples with fixatives/transport medium.
- Replace the cap tightly and invert several times to be sure the tissue floats freely in the liquid and is not stuck to the lid.
- Minimize handling of biopsies; do not use forceps (produces crush artefacts).

Steps:

a. Remove the biopsy core from the needle:

- Use physiologic saline (in a 12 ml syringe fitted with a 25 ga needle)
- GENTLY wash the tissue out of the needle into petridish
- End up with the biopsy core in a small puddle of saline
- Use the saline syringe/needle to wash the needle (with vigor, if needed) before the next pass

Tip: To manipulate the tissue core, use the point of a 25 ga needle and/or the corner of a single-edge razor blade to nudge it around in the saline puddle on the surface of the slide (like tug-boats moving a bigger ship around in a harbor).

- #### b. Promptly assess tissue adequacy (that the cores are of cortex and contain glomeruli) and obtain at least:
- 2 good cores that are ≥ 10 mm long or 3 good cores if they are < 10 mm long.

Divide the total available material for LM, TEM and IF

- Use fresh, clean single-edge razor blade(s) to cut the core(s) transversely
- Separate the pieces for LM, TEM and IF.

Guidelines of how the tissues can be divided for LM, TEM and IF are given in the publication '*Practice guidelines for the renal biopsy*' by Walker and co-authors and The Ad Hoc Committee on Renal Biopsy Guidelines of the Renal Pathology Society, *Modern Pathology*, 2004, 17, 1555–1563.




Diagram to illustrate division of renal biopsy cores in the absence of a dissecting microscope for clinicians using IF. The ends from all cores are taken for TEM with the remainder divided for LM and IF (after Walker).

Tip: Presuming that each portion contains glomeruli, the ‘amount of specimen’ requirements are smallest for TEM, medium for IF, and largest for LM; so, divide the total amount of cortical tissue that is available accordingly.

c. Transfer the pieces into their respective fixative/preservative solutions

- Fix the TEM piece(s) first (should go into the container within 5 minutes)
 - put it (them) in the 3% glutaraldehyde solution
- Fix the LM piece(s) next
 - put it (them) in the clear plastic with 4% formalin
- Fix the IF piece(s) last
 - put it (them) in the Michel’s transport medium

Tip: To transfer the pieces, we usually use a single-edge blade to ‘scoop’ up the piece, which typically clings (due to surface tension) to the ‘wet-with-saline’ edge of the blade, and then we use a pipette containing some fluid from the destination vial to wash the sample off the blade and into the bottle using a new blade and pipette for each tissue transfer task. The last piece(s) on each slide can be washed directly into the destination vial off a corner of the slide.



3. Packing and Return Shipping Instructions

Key aspects:

- Use a commercial courier service to send the kit back to the laboratory.
- Cost of sending the kit back to the laboratory is the responsibility of the clinician/sender.

Please note:

Please be sure the caps are tight on the bottles before placing in the plastic bags.

- **Please contact the laboratory:**
info@mylav.net to confirm expected arrival date of the specimen, name of courier and tracking number. *Please note that the laboratory is closed Sundays and it is therefore best to send the kit back on Monday- Wednesday.*

Address for return shipping:

**Laboratorio MYLAV
Via Sirtori 9,
20017 Rho (Milan)
Italy
tel. +39 02 9308301**



4. Reporting of results and billing:

A preliminary report on the light microscopical findings based upon routine staining techniques (HE, PAS, PAMs and Trichrome) will be send **by e-mail** within 5 days of receipt. A final report which will include descriptions of light, immunofluorescent and electron microscopy and an interpretation will generally be send by e-mail within 2 weeks of receipt.